

powerone)))

Power battery packs
for hearing aids



high level hearing

www.powerone-batteries.com

power one ACCU plus.

Rechargeable NiMH hearing aid batteries.

Your demands

You want hearing aid batteries that you can rely on and that perform well. **power one ACCU plus** are highly efficient batteries based on NiMH technology characterised by the highest quality and reliability.

Your environmental contribution

The quickly recharging batteries are also the ecological first choice. Rechargeable hearing aid batteries are eco-friendly thanks to the usage of environmentally-friendly technologies and packaging materials.

power one batteries are available in size p13 accu, p312 accu and p10 accu.



power one ACCU *plus* meets your expectations:

- › High stability and safe corrosion protection thanks to a stainless steel casing
 - › Additional protection thanks to nanotechnology surface coating
 - › Unlimited functionality under the most varied environmental conditions
 - › Fast-charging capabilities when used with **power one** chargers
 - › No memory effect
 - › Made in Germany
-

Convenient packaging

power one ACCU plus has a packaging which is easy to handle and which safely stores your cells:

- Easy, optimum removal of batteries thanks to a practical rotary dispenser
- Maximum product protection thanks to the quality seal
- The recloseable strap ensures safe storage
- Highly visible colour coding
- Environmentally friendly



power one – fast charger for your mobility.

power one pocketcharger

The **power one** charger can charge one or two **power one ACCU plus** hearing aid batteries.

The **power one pocketcharger** can charge p13 accu, p312 accu or p10 accu within 2.5 hours.

The charger is operated with a rechargeable lithium polymer battery and therefore function independently from a power plug. For longer shelf life of the cells, rechargeable batteries that are still partially charged are automatically discharged before the charging process begins.



power one – Li-Ion.

Rechargeable hearing aid batteries.

The rechargeable hearing aid batteries **power one** Li-Ion are developed individually according to the customer specifications of today and tomorrow. Like all hearing aid batteries from **power one** they are designed for maximum performance and safety.

power one Li-Ion batteries are permanently installed in the hearing aid and cannot be replaced by the user themselves.



Key features:

- Optimised for the limited available space in hearing aids
- Highest possible energy density
- Meets the performance requirements of modern hearing aids and is designed for long and stable use
- Manufactured under strict quality controls
- Made in Germany production



All **power one** batteries are produced in Ellwangen, Germany, in the world's largest and most advanced production facility for hearing aid batteries and are subject to the strictest safety and quality control tests.

Please ensure that the intended charger is used for **power one** brand rechargeable batteries.